SCALPEL HANDLE B.P NO. 5

SKU 1006

[product image not available]

PRODUCT OVERVIEW

Summary	For blades nos. 10, 11, 12 and 15.
Category	General Surgery Instruments, Scal…
Availability	In stock

PRODUCT DESCRIPTION

Uses size 10 - 15 blades Material: Stainless Steel Alias Names: Surgical knife, bard parker handle, B.P handle, dissection knife Surgical Instrument: Knife handles are used to hold various blades to create scalpel. Scalpels are used to make skin incisions or whenever a fine precision cut is necessary. Hospitime Scalpel Handle #5 is used for surgery & anatomical dissection.

PRODUCT TAGS

Scalpel Handle B.P No. 5 | General Surgery Instruments | Scalpel, Blade breaker & Caliper