SCALPEL HANDLE B.P NO. 7

SKU 1007



PRODUCT OVERVIEW

Summary	For blades nos. 10, 11, 12 and 15.
Category	General Surgery Instruments, Scal…
Availability	In stock

PRODUCT DESCRIPTION

Uses size 10 - 15 blades Material: Stainless Steel Alias Names: Surgical knife, BardParker handle, BP handle, Dissection knife Surgical Instrument: Knife handles are used to hold various blades to create scalpel. Scalpels are used to make skin incisions or whenever a fine precision cut is necessary. Hospitime Scalpel Handle #7 is used for surgery & anatomical dissection when precision cutting is needed in a confined space or a deep wound.

PRODUCT TAGS

Scalpel Handle B.P No. 7 | General Surgery Instruments | Scalpel, Blade breaker & Caliper