SKIN HOOK RETRACTOR

SKU 1416

[product image not available]

PRODUCT OVERVIEW

Summary	Sharp or Blunt.
Category	General Surgery Instruments, Retra…
Availability	In stock

PRODUCT DESCRIPTION

Instrument Type: Retractors Prongs: 1 prong Tip Type: Sharp or Blunt Handle: Round Handle Material: Stainless Steel Alias Names: Joseph Skin Hook Surgical Instrument: Hospitime skin hook is a handheld retractor commonly utilized for keeping back the edges of skin during intranasal and pharyngeal methods. The handle is often round and knurled and it tapers to the hook end. This product is straight with a sharp, single prong hook tip. The wide pointedly curved hook towards the jaw is available in a blunt or sharp point. This skin hook is commonly useful for more sensitive surgeries, such as dermatologic procedures. The blunt tip may be better to use on nerves or tendons, and the sharp tip may be ideal for retracting the skin.

PRODUCT TAGS

Skin Hook Retractor | General Surgery Instruments | Retractor