SPENCERWELLS HEMOSTAT FORCEPS

12.5cm Straight | SKU 1178



PRODUCT OVERVIEW

Summary	Straight. Overall length 5in (12.5cm).
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Spencer Wells Instrument Type: Hemostatic Forceps Blade Type: Straight Length: 12.5 cm Tip Type: Serrated Handle: Ratchet & Finger Rings Material: Stainless Steel Alias Names: Spencer-wells artery forceps Surgical Instrument: Hospitime spencer wells hemostatic forceps 12.5cm (straight) are ratcheted hemostatic forceps used to control bloodstream in vessels. They have short, serrated jaws, which can totally seal small blood veins. On larger vessels, they might be utilized to briefly control the blood stream and left to be sutured later.

PRODUCT TAGS

Spencerwells Hemostat Forceps | General Surgery Instruments | Forceps, Clamps & Clips