SPENCERWELLS HEMOSTAT FORCEPS

17.5cm Curved | SKU 1183



PRODUCT OVERVIEW

Summary	Curved. Overall length 7in (17.5cm).
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Spencer Wells Instrument Type: Hemostatic Forceps Blade Type: Curved Length: 17.5 cm Tip Type: Serrated Handle: Ratchet & Finger Rings Material: Stainless Steel Alias Names: Spencer wells artery forceps Surgical Instrument: Hospitime spencer wells hemostatic forceps 17.5cm (curved) are ratcheted hemostatic forceps used to control bloodstream in vessels. They have short, serrated jaws, which can totally seal small blood veins. On larger vessels, they might be utilized to briefly control the blood stream and left to be sutured later.

PRODUCT TAGS

Spencerwells Hemostat Forceps | General Surgery Instruments | Forceps, Clamps & Clips