SPENCERWELLS HEMOSTAT FORCEPS

20cm Curved | SKU 1185



PRODUCT OVERVIEW

Summary	Curved. Overall length 8in (20cm).
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Spencer Wells Instrument Type: Hemostatic Forceps Blade Type: Curved Length: 20 cm Tip Type: Serrated Handle: Ratchet & Finger Rings Material: Stainless Steel Alias Names: Spencer wells artery forceps Surgical Instrument: Hospitime spencer wells haemostatic forceps are ratcheted hemostatic forceps used to control blood stream in vessels. They have short, serrated jaws, which can totally seal small blood veins. On larger vessels they might be utilized to briefly control the blood stream, and left to be sutured later.

PRODUCT TAGS

Spencerwells Hemostat Forceps | General Surgery Instruments | Forceps, Clamps & Clips