
SUTURE CUTTING SCISSORS "HEATH"

15cm | SKU 1075



PRODUCT OVERVIEW

Summary	Suture removal scissors. Serrated j…
Category	General Surgery Instruments, Scisso…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Suture Instrument Type: Scissors Curvature: 'S' Curved Length: 15 cm Tip Type: Serrated & Sharp Ends Handle: Finger Rings Material: Stainless Steel Surgical Instrument: Hospitime suture cutting scissors are essentially utilized for suture removal. Perfect for removing sutures. They should not be utilized for cutting tissues or gauges. they also have serrated jaws.

PRODUCT TAGS

Suture Cutting Scissors "Heath" | General Surgery Instruments | Scissors