
SUTURE TYING DISSECTING FORCEPS

10cm 1X2 Teeth Curved | SKU 1137



PRODUCT OVERVIEW

Summary	1 x 2 teeth. Curved. Overall length 4…
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Corneal Instrument Type: Suture Forceps Curvature: Curved Length: 10 cm Tip Type: 1x2 Teeth & Extra fine with tying platform Handle: Flat with rings Material: Stainless Steel Alias Names: Microdissecting forceps, corneal forceps, Jaffe tying forceps, Mc. Phearson forceps Surgical Instrument: Hospitime suture tying dissecting forceps 10cm (tooth) is an ophthalmic instrument utilized for suturing. The extra fine tips feature a smooth tying platform ideal for use in cataract procedures. The elevated surface holds the suture firmly to help make a knot. This forceps ought not to be used for holding tissues.

PRODUCT TAGS

Suture Tying Dissecting Forceps | General Surgery Instruments | Forceps, Clamps & Clips | OB / GYN Instruments | Forceps & Clamps | Ophthalmic Instruments | Plastic / Dermatology Instruments