
SUTURE TYING DISSECTING FORCEPS

10cm Plain Curved



PRODUCT OVERVIEW

Summary	Plain tying platform. Straight. Overal…
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Corneal Instrument Type: Suture Forceps Curvature: Curved Length: 10 cm Tip Type: Plain & Extra fine with tying platform Handle: Flat with rings Material: Stainless Steel Alias Names: Microdissecting forceps, Corneal forceps, Jaffe tying forceps, Mc. Phearson forceps, Jaffe-Castroviejo suture Forceps Surgical Instrument: Hospitime suture tying dissecting forceps is an ophthalmic instrument utilized for suturing. The extra fine tips feature a smooth tying platform ideal for use in cataract procedures. This forceps ought not to be used for holding tissues.

PRODUCT TAGS

Suture Tying Dissecting Forceps | General Surgery Instruments | Forceps, Clamps & Clips | OB / GYN Instruments | Forceps & Clamps | Ophthalmic Instruments | Plastic / Dermatology Instruments