TISSUE DISSECTING FORCEPS

12.5cm Plain | SKU 1103



PRODUCT OVERVIEW

Summary	Serrated. Overall length 5in (12.5cm…
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Tissue Instrument Type: Tissue Forceps Curvature: Straight Length: 12.5 cm Tip Type: Serrated (Plain) Handle: Thumb Material: Stainless Steel Alias Names: Dissecting Forceps, Thumb forceps, Serrated tissue forceps, Plain dissecting forceps Surgical Instrument: Hospitime tissue dissecting forceps 12.5cm (serrated) are utilized for grasping, holding, and controlling different sorts of tissues in various surgical methods.

PRODUCT TAGS

Tissue Dissecting Forceps | General Surgery Instruments | Forceps, Clamps & Clips