
TISSUE DISSECTING FORCEPS

17.5cm 1X2 Teeth | SKU 1108



PRODUCT OVERVIEW

Summary	1 x 2 teeth. Overall length 7in (17.5…
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Tissue Instrument Type: Tissue Forceps Curvature: Straight Length: 17.5 cm Tip Type: 1x2 Teeth (Tooth) Handle: Thumb Material: Stainless Steel Alias Names: Dissecting Forceps, Thumb forceps, Kocher forceps, Tooth forceps Surgical Instrument: Hospitime tissue dissecting forceps 17.5cm (tooth) are utilized for grasping, holding, and controlling different sorts of tissues in various surgical methods.

PRODUCT TAGS

Tissue Dissecting Forceps | General Surgery Instruments | Forceps, Clamps & Clips