# TISSUE DISSECTING FORCEPS

17.5cm Plain | SKU 1107




## PRODUCT OVERVIEW

| | |
|---|---|
| **Summary** | Serrated. Overall length 7in (17.5cm… |
| **Category** | General Surgery Instruments, Force… |
| **Availability** | In stock |

## PRODUCT DESCRIPTION

Pattern: Tissue Instrument Type: Tissue Forceps Curvature: Straight Length: 17.5 cm Tip Type: Serrated (Plain) Handle: Thumb Material: Stainless Steel Alias Names: Dissecting Forceps, Thumb forceps, Serrated tissue forceps, Plain dissecting forceps Surgical Instrument: Hospitime tissue dissecting forceps 17.5cm (serrated) are utilized for grasping, holding, and controlling different sorts of tissues in various surgical methods.

## PRODUCT TAGS

Tissue Dissecting Forceps | General Surgery Instruments | Forceps, Clamps & Clips