TOWEL CLIP CROSS ACTION

9cm | SKU 1296

[product image not available]

PRODUCT OVERVIEW

Summary	Cross Action. Overall length 3-1/2in…
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Doyen (Cross action) Instrument Type: Towel holding clip Curvature: Curved Length: 9 cm Tip Type: 1x1 sharp prongs Handle: Thumb & cross action Material: Stainless Steel Alias Names: Doyen cross action towel clip Surgical Instrument: Hospitime towel clip is utilized for holding towels or drapes. This item has 1x1 sharp prongs.

PRODUCT TAGS

Towel Clip Cross Action | General Surgery Instruments | Forceps, Clamps & Clips