TOWEL FORCEPS "BACKHAUS"

10cm | SKU 1293

[product image not available]

PRODUCT OVERVIEW

Summary	Overall length 4in (10cm).
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Backhaus Instrument Type: Towel holding forceps Curvature: Curved Length: 10 cm Tip Type: 1x1 sharp prongs Handle: Ratched, Finger Rings Material: Stainless Steel Alias Names: Towel clamp Surgical Instrument: Hospitime Backhaus towel forceps are ratcheted forceps utilized grasping tissues, towels or drapes and holding or reducing small bone fractures. This item is curved and has 1x1 sharp prongs.

PRODUCT TAGS

Towel Forceps "Backhaus" | General Surgery Instruments | Forceps, Clamps & Clips