# TOWEL FORCEPS "BACKHAUS"

12.5cm | SKU 1294

[ product image not available ]

## PRODUCT OVERVIEW

| | |
|---|---|
| **Summary** | Overall length 5in (12.5cm). |
| **Category** | General Surgery Instruments, Force… |
| **Availability** | In stock |

## PRODUCT DESCRIPTION

Pattern: Backhaus Instrument Type: Towel holding forceps Curvature: Curved Length: 12.5 cm Tip Type: 1x1 sharp prongs Handle: Ratched, Finger Rings Material: Stainless Steel Alias Names: Towel clamp Surgical Instrument: Hospitime Backhaus towel forceps are ratcheted forceps utilized grasping tissues, towels or drapes and holding or reducing small bone fractures. This item is curved and has 1x1 sharp prongs.

## PRODUCT TAGS

Towel Forceps "Backhaus" | General Surgery Instruments | Forceps, Clamps & Clips