VOLKMAN RETRACTOR 2 PRONG BLUNT

SKU 1403

[product image not available]

PRODUCT OVERVIEW

Summary	2 Prongs. Blunt. Overall length 8-3/4…
Category	General Surgery Instruments, Retra…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Volkman Length: 22 cm Instrument Type: Retractors Prongs: 2 prongs Tip Type: Blunt Handle: Ring Handle Material: Stainless Steel Alias Names: Volkman ring handle retractor Surgical Instrument: Hospitime volkman retractor is a handheld retractor commonly used for small bone and joint surgeries or procedures. It comes with five different prongs variation (2 prongs to 6 prongs) and in two different tip type (sharp & blunt). Volkman retractor is used to pull wound edges back to allow access.

PRODUCT TAGS

Volkman Retractor 2 Prong Blunt | General Surgery Instruments | Retractor