WAUGH DISSECTING FORCEPS

20cm Plain | SKU 1117



PRODUCT OVERVIEW

Summary	Serrated jaws. Overall length 8in (20…
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Waugh Instrument Type: Tissue Forceps Curvature: Straight Length: 20 cm Tip Type: Fine tip with serrated jaws (Plain) Handle: Thumb Material: Stainless Steel Alias Names: Waugh tissue forceps Surgical Instrument: Hospitime waugh dissecting forceps mostly used in otolaryngology procedures, from ear dissection to dissecting tonsils and placing packing.

PRODUCT TAGS

Waugh Dissecting Forceps | General Surgery Instruments | Forceps, Clamps & Clips | ENT Instruments | Tonsillectomy