
WIRE CUTTING SCISSORS

SKU 1088



PRODUCT OVERVIEW

Summary	Angled jaws. Overall length 5in (12.7…
Category	General Surgery Instruments, Scisso…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Wire Cutting Instrument Type: Scissors Blade Type: Angled Length: 12.7 cm Tip Type: Blunt Ends Handle: Finger Rings Material: Stainless Steel Alias Names: Suture scissor Surgical Instrument: Hospitime wire cutting scissors used for cutting of wires and sutures in preparation for removal. The blades are angled and have serrated jaws.

PRODUCT TAGS

Wire Cutting Scissors | General Surgery Instruments | Scissors